



EUROPEAN COMMISSION
DIRECTORATE-GENERAL
Joint Research Centre

CRLVL13/05XP



Soybean Seeds Sampling and DNA Extraction

Report on the Validation of a DNA Extraction Method from Soybean Seeds

14 May 2007

**Directorate General-Joint Research Centre
Institute for Health and Consumer Protection
Biotechnology & GMOs Unit**

Method development and single laboratory validation:

Bayer CropScience GmbH

Method testing and confirmation:

Community Reference Laboratory for GM Food and Feed (CRL-GMFF)
Biotechnology & GMOs Unit

Drafted by:

E. Grazioli



Report Verification Team:

1) W. Moens




2) M. Querci




Scientific and technical approval:

M. Mazzara



Compliance with CRL Quality System:

S. Cordeil



Authorisation to publish:

G. Van den Eede



Address of contact laboratory:

European Commission, Directorate General-Joint Research Centre

Institute for Health and Consumer Protection (IHCP)

Biotechnology and GMOs Unit – Community Reference Laboratory for GM Food and Feed

Via Fermi 1, 21020 Ispra (VA) - Italy

Content

1. INTRODUCTION.....	4
2. MATERIALS (EQUIPMENT/CHEMICALS/PLASTICWARE).....	4
3. DESCRIPTION OF THE METHODS.....	6
4. TESTING OF THE DNA EXTRACTION METHOD BY THE METHOD DEVELOPER.....	8
5. EXPERIMENTAL TESTING OF THE DNA EXTRACTION METHOD BY THE COMMUNITY REFERENCE LABORATORY FOR GM FOOD AND FEED.....	12
6. CONCLUSION.....	17
7. QUALITY ASSURANCE.....	17
8. REFERENCES.....	17

1. Introduction

This report describes a plant DNA extraction protocol derived from the publicly available “Dellaporta” method ⁽¹⁾. This protocol can be used for the extraction of DNA from soybean seeds and grains ground to powder with a Waring™ blender or with any other appropriate seeds/grains crushing device.

These protocols are recommended to be executed only by skilled laboratory personnel as the procedures comprise the use of hazardous chemicals and materials. It is strongly advised to take particular notice of products safety recommendations and guidelines.

2. Materials (Equipment/Chemicals/Plasticware)

2.1. Equipment

The following equipment is used in the DNA extraction procedure described (equivalents may be used):

1. Waring blender, model 7010S/7010G/7010HS/7010HG or equivalent
2. 70 mm Blender Base (Eberbach Corp. – Cat.No. 8495) for Waring blender or equivalent
3. Micro centrifuge with 18,000 x g for microcentrifuge tubes
4. Table centrifuge (swinging buckets) with 3,000 x g for Falcon tubes
5. Water bath adjustable to 65°C ± 1°C
6. Fluostar Galaxy type 0403 from BMG LabTechnologies
7. PC with Fluostar software (Fluo32)

2.2. Chemicals

The following reagents are used in the DNA extraction procedure described (equivalents may be used):

1. Na₂-EDTA; Titriplex III (Merck Cat. No. 1.08418.1000)
2. Tris-HCl; Tris(hydroxymethyl)aminomethane hydrochloride (USB Cat. No. 22676)
3. NaCl; sodium chloride (Duchefa Cat No. S0520)
4. KAc; potassium acetate (Merck Cat. No. 1.04820.1000)
5. NaAc; sodium acetate (Merck Cat. No. 1.06268.1000)
6. SDS; sodium dodecyl sulphate (BDH Cat. No. 442444H)
7. 2-mercaptoethanol (Sigma Cat. No. M6250)
8. RNase A (Roche Cat. No. 0109-142)
9. Ethanol p.a. (Merck Cat. No. 1.00983.1000)
10. Isopropanol p.a. (Merck Cat. No. 1.09634.2500)
11. Phenol:chloroform:isoamylalcohol (25:24:1) equilibrated at pH 8.0 (Sigma Cat. No. P-2069)

12. Chloroform p.a. (Merck Cat. No. 1.02445.2500)
13. Polyethylene Glycol (MW 8000) (Sigma Cat. No. P2139)

2.3. Solutions

The following buffers and solutions are used in the DNA extraction procedure described:

1. Extraction Buffer pH 8.0

- 100 mM Tris HCl pH 8.0
- 50 mM EDTA
- 500 mM NaCl
- 10 mM 2-mercaptoethanol

2. Tris-EDTA buffer (TE) pH 8.0

- 10 mM Tris HCl pH 8.0
- 1 mM EDTA

3. Tris-EDTA buffer (TE0.1) pH 8.0

- 10 mM Tris HCl pH 8.0
- 0.1 mM EDTA

4. RNase A (10 mg/ml)

5. SDS 20%

6. Ethanol 70%

7. 5M KAc

8. 3M NaAc

2.4. Plasticware

1. 50 ml conical tubes
2. 15 ml conical tubes
3. 2 ml microcentrifuge tubes
4. filter tips

Note: All plasticware has to be sterile and free of DNases, RNases and nucleic acids.

2.5. Precautions

- Phenol, chloroform, isoamyl alcohol, and isopropanol are hazardous chemicals; therefore, all manipulations have to be performed following safety guidelines and under

fume hood.

- It is recommended to use clean containers for Waring blenders for grinding the seed bulk samples.
- All tubes and pipette tips have to be discarded as biological hazardous material

2.6 Abbreviations:

EDTA	ethylenediaminetetraacetic acid
PCR	polymerase chain reaction
RNase A	ribonuclease A
TE	Tris EDTA
Tris	Tris(hydroxymethyl)aminomethane

3. Description of the methods

Sampling:

Sampling approaches for seeds and grains are referred to in technical guidance documents and protocols described in:

- International Organization for Standardizations, Switzerland: ISO standard 6644, ISO standard 13690, ISO standard 5725;
- International Rules for Seed Testing (2004) International Seed Testing Association (ISTA), Switzerland. ISBN 3-906549-38-0;
- USDA-GIPSA (2001) Sampling grains for the detection of Biotech grains <http://www.usda.gov/gipsa/>.

Scope and applicability:

The "Dellaporta-derived" method for DNA extraction described below is suitable for the isolation of genomic DNA from ground soybean seeds and grains. Application of the method to other matrices may require adaptation and possible further specific validation.

Principle:

The basic principle of the DNA extraction consists of first releasing the DNA present in the matrix into aqueous solution and further purifying the DNA from PCR inhibitors. The "Dellaporta-derived" method starts with a lysis step (thermal lysis in the presence of Tris HCl, EDTA, NaCl and 2-mercaptoethanol) followed by isopropanol precipitation and removal of contaminants such as lipophilic molecules and proteins by extraction with phenol:chloroform:isoamylalcohol.

A further DNA precipitate is then generated by using isopropanol, followed by RNase treatment, phenol:chloroform:isoamylalcohol purification and final isopropanol precipitation. The resulting pellet is dissolved in TE buffer.

Seed crushing procedure:

- The seeds/grains are crushed using a Waring blender
- Grind in intervals of 10 seconds for 4 times at maximum speed. Shake between intervals until all powder is loose to improve crushing-procedure
- All seeds/grains should be crushed until a fine powder is obtained. Thorough grinding will also produce a homogenous powder
- Prevent cross-contamination by dust particles between the samples

Soybean seed DNA extraction protocol

1. Transfer 1 g powder into a 50 ml Falcon tube
2. Add 30 ml Extraction Buffer
3. Add 2.1 ml 20% SDS, mix well by inversion
4. Incubate at 65°C for 30 minutes
Note: Mix samples every 10 minutes by inversion
5. Centrifuge for 20 minutes at 3,000 x g
6. Transfer 20 ml supernatant to a new 50 ml Falcon tube using a 25 ml pipette
7. Add 6 ml 5M KAc, shake vigorously for 1 minute
8. Incubate on ice for 30 minutes
Note: Mix samples every 10 minutes by inversion
9. Centrifuge for 20 minutes at 3,000 x g
10. Transfer 20 ml supernatant to a new 50 ml Falcon tube using a 25 ml pipette
11. Add equal volume of isopropanol, mix gently for 1 minute
12. Incubate on ice for 5 minutes
13. Centrifuge for 20 minutes at 3,000 x g
14. Remove supernatant and air-dry the pellet at 37°C until all isopropanol residue is evaporated
15. Dissolve the pellet in 10 ml TE
16. Shake the samples for 1 hour. Make sure the pellet is completely dissolved
17. Add 12 ml phenol:chloroform:isoamylalcohol (25:24:1)
18. Mix well for 1 minute
19. Centrifuge for 20 minutes at 3,000 x g
20. Transfer the upper aqueous phase (8 ml) to a new 50 ml Falcon tube using a 10 ml pipette. Do not disturb the interphase
21. Add 900 µl 3M NaAc
22. Add 6 ml isopropanol
23. Mix gently by inversion for 1 minute
24. Place on ice for 5 minutes
25. Centrifuge for 20 minutes at 3,000 x g to pellet the DNA
26. Remove supernatant and air-dry the pellet at 37°C until all isopropanol residue is evaporated
27. Dissolve the pellet in 1 ml TE

28. Shake the samples for 1 hour. Make sure the pellet is completely dissolved
29. Transfer the DNA solution to a new 2.0 ml microcentrifuge tube
30. Add 10 µl RNase A (10 mg/ml), mix gently and incubate for 20 min at 37°C
31. Add 800 µl phenol:chloroform:isoamylalcohol (25:24:1)
32. Mix well for 1 minute
33. Centrifuge for 10 minutes in a microcentrifuge at maximum speed
34. Transfer the upper aqueous phase (about 900 µl) to a new 2 ml microcentrifuge tube
35. Add 800 µl chloroform
36. Mix well for 1 minute
37. Centrifuge for 10 minutes in a micro centrifuge at maximum speed
38. Transfer the upper aqueous phase (about 800 µl) to a new 2 ml microcentrifuge tube containing 90 µl 3M NaAc
39. Add 600 µl isopropanol
40. Mix gently by inversion for 1 minute
41. Place on ice for 5 minutes
42. Centrifuge for 1 minute in a micro centrifuge at maximum speed to pellet the DNA
43. Remove all supernatant
44. Add 1 ml 70% ethanol to wash the DNA pellet. Make sure the DNA pellet is not stuck to the bottom. Shake the samples for 1 hour
45. Centrifuge for 5 minutes in a micro centrifuge at maximum speed
46. Remove supernatant and air-dry the pellet at 37°C until all ethanol residue is evaporated
47. Add 200 µl TE 0.1 to the DNA pellet
48. Allow pellet to dissolve for about 10 hours at 4°C
49. Shake samples for minimal 3 hours at 4°C
50. Centrifuge for 1 minute in a micro centrifuge at maximum speed
51. Transfer the supernatant to a new microcentrifuge tube

4. Testing of the DNA extraction method by the method developer

4.1 Performance characteristics

A soybean seed bulk sample was ground to fine powder. Six replicated independent DNA extractions were performed on 1 gram flour sub-samples. This experiment was repeated twice under repeatability conditions (within short intervals of time [i.e. days] by the same operator, using the same equipment) resulting in 18 DNA samples.

4.2 DNA concentration, yield and repeatability standard deviation (RSDr)

The concentration of the extracted DNA was determined by fluorescence detection using the PicoGreen dsDNA Quantitation Kit. 1:40 dilutions of each DNA sample were prepared in duplicate and mixed with the PicoGreen reagent. The DNA concentration was determined on

the basis of a nine-point standard curve ranging from 10 ng/ml to 200 ng/ml, using a Fluostar for fluorescence detection. Each concentration was obtained as an average of two readings per sample (Table 1).

Table 1. Results of the 18 (6 X 3) DNA extractions: DNA concentration (ng/μl) and DNA extraction efficiency

DNA extraction	Sample mass (gram)	[DNA] (ng / μl)	Yield (μg DNA / gram sample)
A1	1.0	423.6	84.7
A2	1.0	480.5	96.1
A3	1.0	484.8	97.0
A4	1.0	450.8	90.2
A5	1.0	455.3	91.1
A6	1.0	506.1	101.2
B1	1.0	399.3	79.9
B2	1.0	426.6	85.3
B3	1.0	418.5	83.7
B4	1.0	484.2	96.8
B5	1.0	442.4	88.5
B6	1.0	435.5	87.1
C1	1.0	436.0	87.2
C2	1.0	502.5	100.5
C3	1.0	450.5	90.1
C4	1.0	463.0	92.6
C5	1.0	482.7	96.5
C6	1.0	483.8	96.8

Table 1 lists the DNA concentrations and DNA yield of all 18 DNA extractions. The average DNA concentration was 457 ± 31 ng / μl.


The average DNA extraction efficiency was 91.4 ± 6.1 μg DNA / g of sample, the relative Repeatability Standard Deviation (RSDr) was 6.7%.

Therefore, the soybean DNA extraction procedure is sufficiently efficient and yields adequate amounts of genomic DNA with a sufficient level of repeatability.

4.3 Analysis of DNA fragmentation

Analysis of DNA fragmentation was performed by ethidium bromide-stained agarose gel electrophoresis. One microliter of each undiluted DNA was analyzed on a 1% agarose gel (TBE buffer), including a molecular weight marker (Figure 1).

Figure 1. Agarose gel electrophoresis of eighteen genomic DNA samples extracted from soybean seeds (1 μ l undiluted DNA).



The 18 genomic DNA samples extracted as described above appeared as distinct fluorescent banding patterns migrating through the gel corresponding to high molecular weight DNA. None of the 18 genomic DNA samples showed indications of significant degradation.

4.4 Evidence of the absence of PCR inhibitory compounds

The absence of PCR inhibitory compounds in the DNA preparations was demonstrated by Real-time PCR using the oligonucleotides directed to the endogenous control gene *lectin* on serial dilutions of the DNA preparations.

All DNA extracts were adjusted to a concentration of 40 ng/ μ l. From this sample (named "undiluted sample"), a fourfold serial dilution series was prepared with water (1:4, 1:16, 1:64, 1:256).

To assess the presence of inhibitors, the Ct values (y axis) of the dilution series were plotted against the logarithm of the DNA amount (x axis). By linear regression, a trend line ($y = ax + b$) was calculated, as well as a correlation coefficient, r^2 , as a measure of linearity (Table 2).

The ideal slope value, 'a', (optimal PCR efficiency) then becomes $a = -3.32$ (typically 'a' values between -3.1 and -3.6 indicate excellent PCR efficiencies).

Correlation coefficients of $r^2 > 0.98$ indicate an excellent linear relationship, and thus, equally efficient PCR amplification over the measured dynamic range.

The Ct value for the "undiluted sample" (40 ng/ μ l, 200 ng/reaction) was extrapolated from the equation calculated by linear regression. Subsequently the extrapolated Ct for the undiluted

sample was compared with the measured Ct. Differences between measured and extrapolated Ct of < 0.5 are indicative of low level or absence of inhibition.

The results reported in Table 2 show no evidence of the presence of PCR inhibitors; in fact all 18 Δ Ct values are < 0.5, suggesting absence of PCR inhibitors.

Table 2. Results of the analysis for PCR inhibitory compounds; comparison of extrapolated Ct values versus measured Ct values.

DNA	Slope, a	Intercept, b	Linearity, R ²	Ct _{measured}	Ct _{extrapolated}	Δ Ct (Ct _{measured} - Ct _{extrapolated})
A1	-3.4477	30.1770	0.9997	22.16	22.24	-0.08
A2	-3.5102	30.3226	0.9999	22.35	22.25	0.10
A3	-3.5644	30.2912	0.9995	22.38	22.09	0.29
A4	-3.4918	30.2546	0.9999	22.46	22.22	0.24
A5	-3.4705	30.2414	0.9977	22.34	22.26	0.09
A6	-3.5217	30.4561	0.9999	22.55	22.35	0.20
B1	-3.3889	30.2363	0.9999	22.55	22.44	0.11
B2	-3.4336	30.5378	0.9992	22.79	22.64	0.15
B3	-3.4494	30.6269	0.9999	22.71	22.69	0.02
B4	-3.5080	30.7840	0.9996	23.06	22.71	0.35
B5	-3.3633	30.5680	0.9999	22.93	22.83	0.10
B6	-3.3576	30.5944	0.9995	22.80	22.87	-0.07
C1	-3.4920	30.7847	0.9999	22.79	22.75	0.04
C2	-3.5200	30.7951	0.9996	22.93	22.70	0.24
C3	-3.4440	30.7391	0.9997	22.84	22.81	0.03
C4	-3.4819	30.8506	0.9999	23.08	22.84	0.24
C5	-3.4039	30.6721	0.9999	22.98	22.84	0.15
C6	-3.4822	30.7781	0.9999	22.86	22.77	0.10

5. Experimental testing of the DNA extraction method by the Community Reference Laboratory for GM Food and Feed

The aim of the experimental testing was to verify that the DNA extraction method provides DNA of suitable quantity and quality for the intended purpose.

The DNA extraction method should allow preparation of the analyte in quality and quantity appropriate for the analytical method used to quantify the event-specific analyte versus the reference analyte.

The CRL-GMFF tested the “Dellaporta-derived” method proposed by the applicant on samples of food and feed consisting of ground soybean seeds provided by the applicant.

To assess the suitability of the DNA extraction method for real-time PCR analysis, the extracted DNA was tested using a qualitative PCR run on the real-time PCR equipment.

5.1 Preparation of samples

About 200 g of soybean seed material were ground using a GRINDOMIX GM 200 (Retsch GmbH) mixer.

5.2 DNA extraction

DNA was extracted following the “Dellaporta-derived” method described above; the DNA extraction was carried out on 6 test portions (replicates) and repeated over three different days, giving a total of 18 DNA extractions.

5.3 DNA concentration, yield and repeatability

Concentration of the DNA extracted was determined by fluorescence detection using the PicoGreen dsDNA Quantitation Kit (Molecular Probes).

Each DNA extract was measured twice, and the two values were averaged. DNA concentration was determined on the basis of a five point standard curve ranging from 1 to 500 ng/μl using a Biorad VersaFluor fluorometer.

The DNA concentration for all samples is reported in the Table 3.

Table 3. DNA concentration (ng/ μ l) of eighteen samples extracted in three days: yellow boxes for samples extracted on day 1, green boxes for samples extracted on day 2 and blue boxes for samples extracted on day 3.

Sample	Concentration (ng/ μ l)
1	197
2	188
3	182
4	186
5	190
6	150
1	172
2	147
3	125
4	99
5	91
6	91
1	187
2	159
3	164
4	163
5	177
6	177

✓ DNA concentration (ng/ μ l)

Overall average of all samples:	158 ng/ μ l
Standard deviation of all samples	34.7 ng/ μ l
Coefficient of variation	21.9 %


➤ Yield (total volume of DNA solution: 200 μ l)

Overall average of all samples:	31.6 μ g
Standard deviation	6.9 μ g
Coefficient of variation	21.9 %

5.4 Fragmentation state of DNA

The size of the extracted DNA was evaluated by agarose gel electrophoresis; 8 μ l of the DNA solution were analysed on a 1.0% agarose gel (Figure 1).

Figure 1. Agarose gel electrophoresis of eighteen genomic DNA samples extracted from soybean seeds. Lanes 2-7: samples extracted on day 1; lanes 8-13 samples extracted on day 2; lanes 14-19 samples extracted on day 3; lanes 1 and 20: 1kb DNA Marker.



The eighteen genomic DNA samples extracted as described above appeared as distinct fluorescent banding patterns migrating through the gel corresponding to high molecular weight DNA. None of the DNA samples showed indications of significant degradation ('smearing').

5.4 Purity / Absence of PCR inhibitors

In order to assess the purity and to confirm the absence of PCR inhibitors, the extracted DNA solutions were adjusted to a concentration of 40 ng/μl (hereafter referred as "undiluted" samples).

Subsequently fourfold serial dilutions of each extract were prepared with pure water (1:4, 1:16, 1:64, 1:256) and analysed using a real-time PCR system detecting the target sequence of the endogenous control gene lectin, *Le1*.

The Ct values obtained for "undiluted" and diluted DNA samples are reported in the Table 4.

Table 4. Ct values of undiluted and fourfold serially diluted DNA extracts after amplification of soybean lectin gene, *Le1*. Yellow boxes for samples extracted on day 1, green boxes for samples extracted on day 2 and blue boxes for samples extracted on day 3.

DNA extract	Undiluted (40 ng/μl)	Diluted			
		1:4	1:16	1:64	1:256
1	20.77	22.80	24.90	26.79	29.39
2	22.53	24.63	26.43	28.79	31.11
3	21.07	22.97	25.14	27.03	29.26
4	20.94	22.75	24.78	27.22	29.73
5	21.07	22.99	24.95	26.90	29.38
6	20.94	22.77	24.79	27.02	28.75
1	22.56	24.58	26.86	28.99	31.14
2	20.84	22.99	24.91	27.07	29.31
3	22.68	24.86	26.82	28.87	31.08
4	21.24	23.29	25.33	27.57	29.56
5	21.15	23.13	25.32	27.36	29.67
6	21.15	23.20	25.39	27.53	29.70
1	22.68	24.70	26.78	28.92	31.20
2	22.65	24.88	26.82	28.89	31.38
3	22.78	24.91	26.98	29.22	31.44
4	22.75	24.79	26.78	29.05	31.22
5	22.12	24.15	26.21	28.37	30.65
6	22.53	24.66	26.84	29.06	31.33

Table 5 below reports the comparison of extrapolated Ct values versus measured Ct values for all samples and the values of linearity (R^2) and slope of all measurements.

To measure inhibition, the Ct values of the four diluted samples were plotted against the logarithm of the dilution and the Ct value for the "undiluted" sample (40 ng/μl) was extrapolated from the equation calculated by linear regression.

Subsequently the extrapolated Ct for the "undiluted" sample was compared with the measured Ct. The evaluation is carried out considering that PCR inhibitors are present if the measured Ct value for the "undiluted" sample is suppressed by > 0.5 cycles from the calculated Ct value. In addition, the slope of the curve should be between -3.6 and -3.1.

Table 5. Comparison of extrapolated Ct values versus measured Ct values (amplification of soybean lectin gene, *LeI*)

DNA extraction	R ²	Slope*	Ct extrapolated	mean Ct measured	ΔCt**
1	0.994	-3.6	20.56	20.77	0.21
2	0.996	-3.6	22.28	22.53	0.24
3	0.997	-3.4	20.91	21.07	0.16
4	0.987	-3.9	20.27	20.94	0.67
5	0.995	-3.5	20.77	21.07	0.30
6	0.993	-3.3	20.79	20.94	0.15
1	0.998	-3.6	22.44	22.56	0.12
2	0.997	-3.5	20.79	20.84	0.05
3	0.999	-3.4	22.72	22.68	0.04
4	0.995	-3.5	21.18	21.24	0.06
5	0.999	-3.6	20.96	21.15	0.19
6	0.996	-3.6	21.04	21.15	0.11
1	0.999	-3.6	22.49	22.68	0.20
2	0.996	-3.6	22.60	22.65	0.05
3	1.000	-3.6	22.68	22.78	0.10
4	0.999	-3.6	22.57	22.75	0.19
5	0.999	-3.6	21.93	22.12	0.19
6	1.000	-3.7	22.41	22.53	0.12

Note: In yellow boxes samples from 1 to 6 extracted on day 1; in green boxes samples from 1-6 extracted on day 2; in blue boxes samples from 1-6 extracted on day 3.

*The expected slope for a PCR with 100% efficiency is -3.32

**delta Ct = abs (Ct extrapolated - Ct measured)

All ΔCt values of extrapolated versus measured Ct are < 0.5, except the sample number 4 extracted on day 1, with a value of 0.67.

R² of linear regression is > 0.99 for all DNA samples except one slightly below (0.987).

Curve slopes are all between -3.1 and -3.6, with two exceptions: sample 4 extracted on day 1, with a value of -3.9 (same sample showing a ΔCt of 0.67) and sample 6 extracted on day 3, with a value of -3.7.

6. Conclusion

The data reported confirm that the extraction method, applied to soybean seeds provided by the applicant, produces DNA of suitable quantity and quality for subsequent PCR based detection applications. The method is consequently applicable to samples of soybean seeds provided as samples of food and feed in accordance with the requirements of Annex I-2.C.2 to Commission Regulation (EC) No 641/2004.

7. Quality assurance

The CRL-GMFF carries out all operations according to ISO 9001:2000 (certificate number: CH-32232) and ISO 17025:2005 (certificate number: DAC-PL-0459-06-00) [DNA extraction, qualitative and quantitative PCR in the area of Biology (DNA extraction and PCR method validation for the detection and identification of GMOs in food and feed materials)]

8. References

1. Dellaporta *et al.* [1983] Plant Molecular Biology Reporter 1, nr 4, pp. 19-21.
2. Sambrook J. and D. W. Russell. Molecular Cloning. A laboratory manual. Third edition. 2001. Cold Spring Harbor Laboratory Press.

